

# RNase Inhibitor V2

## Enhance the Quality and Reliability of Your Single-Nuclei Sequencing Experiments

### Description

Gene expression assays that utilize single nuclei as input require inhibition of endogenous RNase activity during the nuclei isolation process. RNase Inhibitor V2 has been optimized to inhibit endogenous RNase activity, including RNases A, B, and C, during nuclei isolation and purification from solid tissue, thereby ensuring the integrity and stability of RNA within single nuclei for gene expression detection. The RNase Inhibitor V2 is designed to be used with the Singulator platform, with packaging sizes matching the number of isolations and subsequent processing steps. Our V2 RNase Inhibitor exhibits improved performance over our V1 RNase Inhibitor, ensuring the success of your single-nuclei genomics experiments. The RNase Inhibitor V2 is compatible with several downstream applications, including the 10x Genomics' Chromium Gene Expression and Multiome Assays and Parse Evercode™ WT Assays.


### Product Features

- Enhanced Activity:** RNase Inhibitor V2 demonstrates enhanced activity over its predecessor, V1.
- Robust RNase Inhibition:** Robust protection against RNases A, B, and C during nuclei isolation, ensuring RNA integrity with its high stability and maximum inhibition.
- Seamless Integration with the Singulator Platform:** RNase Inhibitor V2 performs with Singulator Nuclei Isolation Kits, sold in sizes that cover both isolation and post-processing needs.
- Economically Priced:** RNase Inhibitor V2 makes automated single-nuclei isolation affordable without sacrificing quality.


### Performance

The RNase Inhibitor V2 was validated using fresh and frozen samples into single-nuclei suspensions suitable for downstream genomic analysis.

The RNase Inhibitor V2 was compared to the RNase Inhibitor V1 in nuclei isolated from a pancreatic ductal adenocarcinoma (PDAC) mouse model tumor, a tissue with a high level of endogenous RNase activity. Nuclei were isolated using the Standard Nuclei Isolation Protocol on the Singulator Platform in combination with the NIC+ cartridges. Gene expression observed using the 10x Genomics Chromium Next GEM Single Cell 3' Kit v3.1. Performance was assessed by measuring the unique gene number and total UMI counts of samples isolated with RNase Inhibitor V2 compared to V1 (**FIGURE 1**). Samples prepared with RNase Inhibitor V2 had 15.6K genes detected compared to 14.7K with V1, and total UMIs detected with V2 were 14.1M versus 11.6M with V1. This indicates V2 prepared samples have more unique genes with higher counts compared to V1 prepared samples.



**FIGURE 1: Gene Expression difference between mouse PDAC nuclei isolated with V1 and V2 RNase Inhibitor.** This graph illustrates that for all threshold values (K), the number of genes with counts greater than K is consistently higher in V2 compared to V1.



## Specifications

Concentration: 40U/uL  
Unit definition: One unit of RNase inhibitor is defined as the amount of RNase Inhibitor required to inhibit activity of 0.375 ng of RNase A by  $\geq 95\%$

ORDERING INFORMATION:	
Item Description	Part Number
RNase Inhibitor V2, 8 samples	100-288-916
RNase Inhibitor V2, 24 samples	100-291-086
RNase Inhibitor V2, 96 samples	100-291-195
Nuclei Isolation Bundle with RNase Inhibitor V2 (24 samples)	100-288-807
NIC+ Isolation Bundle with RNase Inhibitor V2 (24 samples)	100-289-152
FFPE Nuclei Isolation Bundle and RNase Inhibitor (24 samples)	100-284-221

## Contact Us

→ For general information and technical resources:



<https://s2genomics.com>

→ For sales information or to purchase products:



[inquires@s2genomics.com](mailto:inquires@s2genomics.com)

### Phone

1-925-292-8243

### Address

S2 Genomics, Inc  
7683 Southfront Rd. Suite 200  
Livermore, CA 94551  
USA

### LinkedIn

[Linkedin.com/company/s2genomics](https://www.linkedin.com/company/s2genomics)

### Twitter

[Twitter.com/s2genomics](https://twitter.com/s2genomics)

THIS PRODUCT IS FOR RESEARCH USE ONLY.