



## Singulator™ 100 Nuclei Protocol Submission

Please fill out form as thoroughly as possible. For additional questions and support, email [community@s2genomics.com](mailto:community@s2genomics.com). For more information on the Singulator™ 100 and single-cell processing, head over to [www.s2genomics.com](http://www.s2genomics.com).

### General Information | Study Identification

Protocol Name:

Investigator Name:

Investigator Email:

Secondary Investigator(s):

Secondary Email(s):

Tissue Species:

Tissue Type:

Tissue State:

Mass (mg):

Pre-Singulator™ 100 Processing | Run Summary:

### Singulator™ 100 Nuclei Protocol Parameters

Reagents:

Custom Formulation:

Protocol Type:

Auto Mince:    Yes    No

Incubation Time:

Incubation Temperature:

Mixing Type:

Mixing Speed:

Disruption Type:

Disruption Speed:

### Post-Singulator™ 100 Processing

Centrifuge Time & Speed:

Additional Cleanups/Notes: